

WE COLLECT THE CELLS YOU SEARCH FOR

GILUPI
MEDICINE. PERSONALIZED.

GILUPI GmbH Hermannswerder 20a 14473 Potsdam Germany

Attention:

If you are using **BD Venflon™ Pro Cannula** for the application of the GILUPI CellCollector®

Affected Products:

BD Venflon™ Pro IV Cannula

BD Venflon Pro 20GA 1,1 mm x 32 mm, Order 393204

BD Venflon™ Pro Safety (VPS) Needle Protected IV Cannula

BD Venflon Pro Safety 20GA 1,1 mm x 32 mm PZN 13156910; Order 393224

BD Venflon Pro Safety 20GA 1,1 mm x 32mm, INSTAFLASH, Order 393281


Since May 2021 there are Recalls ongoing!

Statement from BD:

“Based on the current situational assessment, BD has elected to initiate a global recall affecting all ethylene oxide (EtO) sterilized BD Venflon™ Pro IV Cannula. This product recall does not impact BD Venflon™ Pro IV Cannula which has been sterilized by radiation (E-beam).

Note: The products impacted by this recall are not sold in North America; this recall does not impact the US/Canada markets.”

“Please immediately review your inventory to determine if you have potentially affected devices. Follow the links to the searchable database below to search by product lot number(s) or product catalog number(s).”



<http://www.bd.com/MDS-21-4111>



<https://www.bd.com/MDS-20-3801>

Please do not use the GILUPI CellCollector® with the affected BD products!

Written by PRRC on May 7, 2021 on behalf of GILUPI GmbH